**Descrición:**

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**Description:**

**Alpha-Flow**<sup>®</sup> is a radiopaque, fluoride containing, light cure, resin-based micro-hybrid flowable composite. The organic matrix is based on Bis-GMA resin and contains a 38% by volume blend of inorganic fillers with particle range from 0.02 - 0.7 microns. This special combination provides a material with extreme versatility for use in a wide range of procedures, including Class III and V cavity figurations and small Class I and II restorations. This material is not suited for restoration of occlusal surfaces. **Alpha-Flow**<sup>®</sup> is available in our Standard Kit (contains shades A2, A3, B2, and Universal Opaquer) or Custom Kit (available in most popular Vita<sup>®</sup> shades). **Alpha-Flow**<sup>®</sup> now contains fluorescence to provide a duplication of the natural tooth substance.

**Indications for Use:**

- Small Class I and II restorations (non-occlusal surface)
- Class III and V restorations
- Block out of small undercuts
- Fill tunnel preps
- Marginal defect repairs
- Repair composite and ceramic veneers
- Re-conform temporary restorations

**Contraindication:**

**Alpha-Flow**<sup>®</sup> is contraindicated in patients with a history of severe allergic reaction to methacrylate resins.

**Warnings:**

- Avoid contact of etchant with soft tissue, eyes and skin. **Alpha-Etch37**<sup>®</sup> contains phosphoric acid, which is corrosive. In case of contact with soft tissue, eyes or skin immediately flush with water for at least 15 minutes. Seek immediate medical attention.
- Do not use in the presence of eugenol based materials as complete setting may not occur.
- Applicator tips are disposable and intended for single patient use.

**General Precautions:**

- U. S. Federal Law restricts this device to sale by or on the order of a dentist.
- Review instructions and SDS (www.dentaltech.com/resource) prior to beginning treatment. Use only as directed.

**Adverse Reaction:**

Allergic or pre-existing skin, eye and respiratory conditions may be aggravated by exposure to this product.


**Instructions for Use:**

1. Prepare the tooth. Isolate and dry the tooth thoroughly.
2. If using **Alpha-Etch37**<sup>®</sup> etch enamel for 15 seconds or if applied to dentin etch for 10 seconds; never leave etchant on tooth for more than 60 seconds. Lightly agitating the etchant during its application can enhance its action.
3. Rinse thoroughly with water.
4. If using **Alpha-Bond**<sup>®</sup> Light Cure Bonding Resin, thoroughly dry the etched tooth. Dispense two drops of **Alpha-Bond**<sup>®</sup> Light Cure Bonding Resin onto a dispensing pad or other suitable surface. Thoroughly saturate a brush tip and generously apply one uniform layer onto the prepared tooth. Light cure the bonding resin by applying a dental halogen light for 15 - 20 seconds. The surface should appear glossy. If not, repeat application. If using a light cure bonding adhesive from another manufacturer, follow their instructions to determine if the surface of the tooth is to be thoroughly dried or remain moist.
5. Select the appropriate shade of **Alpha-Flow**<sup>®</sup>. Apply **Alpha-Flow**<sup>®</sup> directly onto the prepared tooth in increments no larger than 2 mm thickness.
6. Using a visible light-curing unit, with a minimum output of 600 mW/cm<sup>2</sup>, cure composite material for at least 40 seconds. You can expect the depth of cure for light shades to be 5.0 mm and for dark shades 4.0 mm.
7. Finish and polish using standard techniques. A composite surface sealant can also be used to improve surface smoothness and seal micro cracks and surface imperfections caused by the finishing procedures.

**Storage:**

- **Alpha-Flow**<sup>®</sup> is formulated to be used at room temperature. See outer package for expiry date. Do not freeze. Recommended storage conditions are between 2°C/36°F to 27°C/80°F, bringing to room temperature before using (23°C/73°F).
- Do not expose restorative materials to elevated temperatures or intense light.
- Always replace the syringe cap after each use.

\* Vita is a registered trademark of VITA ZAHNFABRIK H. Rauter GmbH & Co. KG Germany.

|    |  |
|----|--|
| DE | LICHTHÄRTEND HARZBASIERTES FLIESSFÄHIGES MIKRO-HYBRID-KOMPOSIT MIT FLUORID             |
| ES | FOTOCURADO COMPUESTO FLUIDO MICROHIBRIDO A BASE DE RESINA CON FLUORURO                 |
| FR | PHOTOPOLYMÉRISATION COMPOSITE MICROHYBRIDE FLUIDIFIABLE À BASE DE RÉSINE AVEC FLUORURE |
| IT | FOTOPOLIMERIZZAZIONE COMPOSITO FLUIDO MICRO-IBRIDO A BASE DI RESINA CON FLUORURO       |
| LT | KIETINIMAS ŠVIESA MIKROHIBRIDINIS TAKUSIS DERVOS PAGRINDO KOMPOZITAS SUFLUORIDU        |
| PL | UTWARDZANY ŚWIATŁEM PŁYNNY KOMPÓZYT MIKROHYBRYDOWY NABAZIE ŻYWICY ZAWIERAJĄCY FLUOREK  |
| RO | FOTOPOLIMERIZABIL MATERIAL COMPOZIT FLUID MICROHIBRID PE BAZĂ DE RĂȘINĂ CU FLUOR       |
| TR | FLORÜRLÜ REZİN BAZLI MIKROHİBRİD AKIŞKAN KOMPOZİT                                      |

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\* Vita, VITA ZAHNFABRIK H. Rauter GmbH & Co. KG, Almanya firmasının tescilli ticari markasıdır.