Descripción:

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Description:

Alpha-dent[®] Zinc Phosphate Cement is a powder/liquid system for manual mixing. **Alpha-dent**[®] Zinc Phosphate Cement has low thermo conductivity making it an excellent insulator under permanent restorations. The powder is available in two shades; light yellow or white.

Indications for Use:

- Cementation of inlays, onlays, crowns, bridges, and fixed orthodontics.
- Intermediate base material under amalgam or composite.
- Temporary restorations

Contraindications:

Alpha-dent[®] Zinc Phosphate Cement is contraindicated in direct and indirect pulp capping.

Warnings:

- Avoid contact of **Alpha-dent**[®] Zinc Phosphate Cement with soft tissue, eyes and skin. It contains phosphoric acid which is corrosive. In case of contact with soft tissue, eyes or skin immediately flush with water for at least 15 minutes. Seek immediate medical attention.

General Precautions:

- U.S. Federal Law restricts this device to sale by or on the order of dentist.
- Review instructions and SDS (www.dentaltech.com/resource) prior to beginning treatment. Use only as directed.

Adverse Reaction:

Allergic or pre-existing skin, eye and respirator conditions may be aggravated by exposure to this product.

Instructions for Use:

1. Temperature influences working and setting times. At higher temperatures, shorter working time and setting time will result. If a longer working and setting time is desired, mix the cement on a cold, dry glass slab. Any condensation that occurs on the slab will not affect the quality of the cement. The temperature range for the preparation should be 18-24°C (65-75°F). Adding more powder will produce a thicker mix to be used as an intermediate base and temporary restorations . Additional powder will decrease working and setting times, however, once mixing has begun do not add liquid.
 2. Prepare the tooth. Isolate and dry the tooth thoroughly.
 3. To obtain the desired consistency for cementation, dispense the powder and liquid next to each other on a glass slab, or other suitable surface, using the following proportions: 3 drops (0.13 grams) of **Alpha-dent**[®] Zinc Phosphate Liquid to 1 level scoop of **Alpha-dent**[®] Zinc Phosphate Powder (0.20 grams).
 4. Divide the powder into four-equal parts or quarters. Divide one of these quarter portions into two parts, or two eighths. Then divide one of these eighth portions into two parts, or two sixteenths. You should now have 6 portions (1/4, 1/4, 1/4, 1/8, 1/16, 1/16).
 5. Using a plastic spatula, mix one of the 1/16 portions of powder into the entire liquid for 10 seconds. Add the other 1/16 and mix for 10 seconds. Add the 1/8 portion and mix for 10 seconds. Then add consecutively two of the quarter portions mixing for 15 seconds after each addition. Add the final quarter portion and mix the entire amount for 30 seconds. This will give a total mixing time of 90 seconds. If cement is too viscous it should be discarded, and the mixing procedure restarted.
 6. Apply **Alpha-dent**[®] Zinc Phosphate Cement in a uniform layer. The working time is approximately 2 minutes from the start of mixing so placement of the cement should be done immediately. Complete setting of **Alpha-dent**[®] Zinc Phosphate Cement will occur 5 minutes and 30 seconds from the end of mixing (net setting time).
 7. Excess cement should be removed at once with cotton or other suitable instrument or when the cement has completely hardened.
- Storage:**
- Store material between 2°C/36°F to 27°C/80°F, bringing to room temperature before using (23°C/73°F). Do not freeze. See outer package for expiration date.
 - Keep containers closed when not in use.

DE	ZINKPHOSPHATZEMENT	
ES	CEMENTO DE FOSFATO DE ZINC	
FR	CIMENT AU PHOSPHATE DE ZINC	
IT	CEMENTO AL FOSFATO DI ZINCO	
LT	CINKO FOSFATO CEMENTAS	
NL	ZINKFOSFAATCEMENT	
PL	CEMENT FOSFORANU CYNKU	
RO	CIMENT FOSFAT DE ZINC	
TR	ÇİNKO FOSFAT SİMAN	
		PM-2586 Rev011

6901 N. HAMLIN AVENUE LINCOLNWOOD, ILLINOIS 60712 USA www.dentaltech.com	0459 EMERGO EUROPE Prinsessegracht 20 2514AP, The Hague The Netherlands

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