**Descripción:**

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**Description:**

Alpha-Glaze® surface sealant is a light cure fluoride containing unfilled methacrylate-based resin with low viscosity used as a finishing sealant for composite and resin cements. It is formulated to penetrate and seal the micro-cracks and surface imperfections that are created by finishing procedures.

Indications for Use:

- Seals composite restoration surfaces.
- Increases surface hardness.
- Increases the resistance to wear.
- Seals tooth surface that helps protect the enamel and makes paste clean up easier prior to applying Orthodontic Appliances.

Contraindications:

Alpha-Glaze® is contraindicated in patients with a history of severe allergic reaction to methacrylate resins.

Warning:

- Avoid contact of **Alpha-dent**® etchant and **Alpha-Etch37**® with soft tissue, eyes and skin. Both contain phosphoric acid, which is corrosive. In case of contact with soft tissue, eyes or skin immediately flush with water for at least 15 minutes. Seek immediate medical attention.

General Precautions:

- U. S. Federal Law restricts this device to sale by or on the order of a dentist.
- Review instructions and SDS (www.dentaltech.com/resource) prior to beginning treatment. Use only as directed.
- Always replace cap after each use.

Adverse Reaction:

Allergic or pre-existing skin, eye and respiratory conditions may be aggravated by exposure to this product.

Instructions for Use:**Finishing Sealant for Composite and Resin Cements**

1. Following the final adjustment, contouring, finishing and polishing of the composite restoration, rinse the tooth completely to remove all debris
2. Apply **Alpha-dent**® etchant or **Alpha-Etch37**® to the prepared restoration surface, and approximately 1-2mm of the surrounding enamel. For full coverage restorations, etch the entire surface of the restoration. Etch the surface for 15 seconds. Lightly agitating the etchant during its application can enhance its action. Never leave the etchant on the tooth for more than 60 seconds.
3. Rinse thoroughly with water, dry the tooth and restoration with oil and moisture-free air. The dry etched surface should appear dull and frosty white. If the surface that was etched does not, then repeat step 2.
4. Place small amount of **Alpha-Glaze**® in a mixing well. Using a clean brush tip, apply a thin layer to the previously etched tooth and restorative surfaces.
5. Using a light-curing unit (minimum output of 600mW/cm2), light cure **Alpha-Glaze**® for 20 seconds. Check occlusion. Typically, the occlusion will not require adjustment due to the minimal film thickness of **Alpha-Glaze**®.

Finishing Sealant for Orthodontics

1. Prophylaxis the tooth structure to be bonded, rinse thoroughly and dry the teeth.
2. Apply **Alpha-dent**® etchant or **Alpha-Etch37**® to the prepared enamel for 15 seconds. Lightly agitating the etchant during its application can enhance its action. Never leave the etchant on the tooth for more than 60 seconds.
3. Rinse thoroughly with water, and dry the tooth with oil and moisture-free air. The dry etched surface should appear dull and frosty white. If the surface that was etched does not, then repeat step 2.
4. If using **Alpha-dent**® Bracket Bonding Resin, it should be applied at this step, following the instructions before proceeding.
5. Using a brush or other suitable instrument, apply a thin layer of **Alpha-Glaze**® onto the dry labial surface.
6. Using a light-curing unit (minimum output of 600mW/cm2), Light cure **Alpha-Glaze**® for 20 seconds.
7. Tooth is ready for bonding of appliance.

Storage:

- Recommended storage conditions are between 2°C/36°F to 27°C/80°F, bringing to room temperature before using (23°C/73°F). Do not freeze. See outer package for expiration date.
- Do not expose to elevated temperatures or intense light.
- Keep container closed when not in use.

	KOMPOSITE OBERFLÄCHENVERSIEGLER LICHTHÄRTEND
	COMPOSITE SELLADOR DE SUPERFICIE FOTOCURADO
	RÉSINE COMPOSITE PRODUIT DE SCELLEMENT DE SURFACE PHOTOPOLYMÉRISATION
	COMPOSITO SIGILLANTE SUPERFICIALE FOTOPOLIMERIZZAZIONE
	KOMPOZITAS PAVIRŠIAUS SANDARIKLIS KIETINAMAS ŠVIESA
	COMPOSIT PPPERVLAKTE-SEALANT LICHTUITHARDEND
	KOMPOZYT USZCZELNIACZ DO POWIERZCHNI UTWARDZANY ŚWIATŁEM
	MATERIAL COMPOZIT MATERIAL DE SIGILARE PENTRU SUPRAFEȚE FOTOPOLIMERIZABIL
	KOMPOZİT YÜZEY ÖRTÜCÜ İŞIKLI POLİMERİZASYON
	2°C/36°F 27°C/80°F

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