ES

Descripción:

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EN

Description:

Alpha-dent[®] Polycarboxylate Cement is specially formulated for superior strength, resistance to dissolution, and oral retention. Other excellent physical characteristics exhibited by **Alpha-dent**[®] Polycarboxylate Cement include; low film thickness for complete seating of appliances, ionic bonding to enamel and dentin, and mild to pulp and tissue.

Indications for use:

- Cementation of inlays, onlays, crowns, bridges, and fixed orthodontic bands.
- Cavity lining under all restorative materials or following root canal treatment.
- Temporary fillings in deciduous teeth.

Contraindications:

Alpha-dent[®] Polycarboxylate Cement is contraindicated in patients with a history of dermatitis and respiratory medical conditions.

Warning:

- Avoid contact of **Alpha-dent**[®] Polycarboxylate Cement with soft tissue, eyes and skin. It contains polyacrylic acid which is corrosive. In case of contact with soft tissue, eyes or skin immediately flush with water for at least 15 minutes. Seek immediate medical attention.

General Precautions:

- U.S. Federal Law restricts this device to sale by or on the order of a dentist.
- Review instructions and SDS (www.dentaltech.com/resource) prior to beginning treatment. Use only as directed.

Adverse Reaction:

Allergic or pre-existing skin, eye and respiratory conditions may be aggravated by exposure to this product.

Instructions for Use:

The recommended temperature range for preparation is 18°C/64°F to 24°C/75°F.

Dispensing Material

1. Tap powder bottle on a firm surface to condense powder.
2. Scoop powder with the measuring scoop. Level off excess powder with a spatula, leaving a flat surface.
3. Place measured powder onto the mixing pad or a glass slab by tapping the end of the measuring scoop with a spatula.
4. When dispensing the liquid, hold the dropper bottle at a slight angle to vertical and dispense the number of drops required. Squeeze bottle to produce the drop, and release pressure after each drop is dispensed. Remove excess liquid from the dropper tip with a damp cloth to avoid clogging of the nozzle.

Dispensing Ratio

1. For cementation: 1 measured scoop of powder (0.22g) to 2 or 3 drops of liquid (0.09g - 0.14g), depending on the consistency preferred.
2. For linings and fillings: 1 measured scoop of powder (0.22g) to 1 drop of liquid (0.045g).

Mixing

Using a spatula or other suitable instrument, incorporate all of the powder rapidly into the liquid. Mix until a creamy consistency is reached for cementation material and a putty-like consistency for lining/filling material. The mixing time is 30 seconds. The working time is 1 minute and 45 seconds from the start of mixing.

Application

Apply the cement immediately upon completion of the mix, while the appearance is glossy. Do not use after it starts to "cobweb". The setting time is 7 minutes maximum at 23°C/73°F. The net setting at 37°C/99°F is 3 minutes from the end of mixing.

Finishing

Do not disturb the material when it is in the intermediate rubbery stage. Excess material may be removed at once or when the cement has completely hardened.

Clean-up

Due to the adhesive nature of **Alpha-dent**[®] Polycarboxylate Cement; spatulas, glass slabs, and placement instruments should be wiped clean under running water immediately after use to ease clean-up.

Storage:

- Store material between 2°C/36°F to 27°C/80°F, bringing to room temperature before using (23°C/73°F). Do not freeze. See outer package for expiration date.
- Keep containers closed when not in use.

**DENTAL
TECHNOLOGIES**alpha-dent[®]**POLYCARBOXYLATE
CEMENT**

DE POLYCARBOXYLATZEMENT

ES CEMENTO DE POLICARBOXLIATOS

FR CIMENT POLYCARBOXYLATE

IT CEMENTO AL POLICARBOSSILATO

LT POLIKARBOKSILATO CEMENTAS

NL POLYCARBOXYLAATCEMENT

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RO CIMENT POLICARBOXILAT

TR POLİKARBOKSİLAT SİMAN

**CORROSIVE**

PM-2943 Rev009

**DENTAL
TECHNOLOGIES**6901 N. HAMLIN AVENUE
LINCOLNWOOD, ILLINOIS 60712 USA
www.dentaltech.com

0459

EMERGO EUROPE

Prinsessegracht 20

2514 AP, The Hague

The Netherlands