ES

**Descripción:**

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EN

**Description:**

**Alpha-dent**® Glass Ionomer Cement is a specially formulated fluoride releasing cement which provides: superior strength, resistance to dissolution, and oral retention. It adheres to dentin, enamel, plastic and non-precious metals.

**Indications for Use:**

- For permanent cementation of non-precious metal inlays, crowns, bridges, and orthodontic bands.
- Intermediate base material under amalgam or composite.

**Contraindications:**

**Alpha-dent**® Glass Ionomer Cement is contraindicated in direct and indirect pulp capping.

**Warnings:**

**Alpha-dent**® Glass Ionomer Cement contains polyacrylic acid, which is corrosive. Avoid contact with eyes. In case of contact with soft tissue, eyes or skin immediately flush with water for at least 15 minutes. Seek immediate medical attention.

**General Precautions:**

- U.S. Federal Law restricts this device to sale by or on the order of a dentist.
- Review instructions and SDS (www.dentaltech.com/resource) prior to beginning treatment. Use only as directed.
- Excess **Alpha-dent**® Glass Ionomer Cement should be removed with a cotton roll immediately before it sets. All surrounding areas such as gingiva and adjacent teeth can be lubricated with petroleum jelly to prevent adhesion. These areas should be cleaned before cement reaches its final set.

**Allergic Reaction:**

Allergic or pre-existing skin, eye and respiratory conditions may be aggravated by exposure to this product.

**Instructions for Use:**

- Thoroughly clean teeth to be bonded with **Alpha-Pro**® No Added Color/Flavor/Fluoride Prophyl Paste. Prophyl paste which contains oil or fluoride should not be used as it may leave a residue which can inhibit bonding.
- Rinse teeth thoroughly with large amount of water and dry with oil/moisture-free compressed air. Provide suitable moisture control to isolate area and to prevent moisture contamination during initial setting.
- Prior to measuring, shake powder bottle to uniformly blend the powder.
- Dispense 1 full scoop of powder and 2 drops of liquid on a mixing pad or glass slab. When dispensing the liquid, hold the dropper bottle in a slight angle to vertical and dispense the number of drops required. Squeeze bottle to produce the drop, and release pressure after each drop is dispensed. Remove excess liquid from the dropper tip with a damp cloth to avoid clogging of the nozzle.
- Incorporate the powder into the liquid in small increments. Total mixing time should not exceed 30 seconds.
- Apply **Alpha-dent**® Glass Ionomer Cement in a uniform layer. Placement of mixed cement should be done immediately while mix is still glossy. If cement becomes dull, it should be discarded and a new mix prepared.
- The working time of **Alpha-dent**® Glass Ionomer Cement is approximately 3 minutes at 23°C/73°F. If a longer working time is needed, mix the cement on a cool, dry glass slab.
- Keep cement free from saliva or water for 5 minutes. At 37°C/99°F, complete setting of **Alpha-dent**® Glass Ionomer Cement will occur 5 minutes and 30 seconds from the end of mixing (net setting time).

**Storage:**


- Store material between 2°C/36°F to 27°C/80°F, bringing to room temperature before using (23°C/73°F). Do not freeze. See outer package for expiration date.
- Keep containers closed when not in use.

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