#### Descripción:

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#### Description:

**Alpha II® AP** is a light cure, resin-based, micro-hybrid composite with Fluoride. **Alpha II® AP** is available in six Vita™ shades plus three opaquers and an incisal shade. The organic matrix oligomer is based on Bis-GMA resin. **Alpha II® AP** contains a 53% by volume blend of inorganic fillers with particle range from 0.02-2 microns. **Alpha II® AP** is a radio-opaque restorative material that allows easy identification.

#### Indications for Use:

**Alpha II® AP** is recommended for use in Class I and Class II restorations and is aesthetically accurate for use in Class III, IV, and V restorations.

#### Contraindications:

**Alpha II® AP** is contraindicated in patients with a history of severe allergic reaction to methacrylate resins.

#### Warnings:

- Avoid contact of etchant with soft tissue, eyes, and skin. **Alpha-Etch37®** contains phosphoric acid, which is corrosive. In case of contact with soft tissue, eyes or skin immediately flush with water for at least 15 minutes. Seek medical attention.
- Do not use in the presence of eugenol based materials as complete setting may not occur.

#### General Precautions:

- U. S. Federal Law restricts this device to sale by or on the order of a dentist.
- Review instructions and SDS (www.dentaltech.com/resource) prior to beginning treatment. Use only as directed.

#### Adverse Reactions:

Allergic or pre-existing skin, eye and respiratory conditions may be aggravated by exposure to this product.

#### Instructions for use:

- Prepare the tooth. Isolate and dry the tooth thoroughly.
- If using **Alpha-Etch37®** etch enamel for 15 seconds or if applied to dentin etch for 10 seconds; never leave etchant on tooth for more than 60 seconds. Lightly agitating the etchant during its application can enhance its action.
- Rinse thoroughly with water.
- If using **Alpha-Bond®** Light Cure Bonding Adhesive, the surface must remain moist. Dispense two drops of **Alpha-Bond®** Light Cure Bonding Adhesive onto a suitable surface. Thoroughly saturate an applicator tip and generously apply two coats onto the tooth, remove excess adhesive by suctioning being careful to keep a continuous coating on the tooth. This should be followed by gently applying air for at least 5 seconds to evaporate the solvent, making sure you do not apply too much air to thin the adhesive. When properly dried the adhesive layer will not move with air and the surface should appear shiny. If the surface is not shiny it may not be properly sealed and additional applications are recommended. Light cure the adhesive for 10 seconds (additional curing time may be required if output of curing unit is less than 600 mW/cm<sup>2</sup>).
- Select the appropriate shade of **Alpha II® AP**. Apply **Alpha II® AP** directly onto the prepared tooth in increments no larger than 2mm thickness.
- Using a visible light-curing unit, with a minimum output of 600 mW/cm<sup>2</sup>, cure each incremental layer for at least 40 seconds.
- Finish and polish using standard techniques. A composite surface sealant can also be used to improve surface smoothness and seal micro cracks and surface imperfections caused by the finishing procedures.


#### Storage:

- Alpha II® AP** is formulated to be used at room temperature. See outer package for expiry date. Do not freeze. Recommended storage conditions are between 2 °C/36 °F to 27 °C/80 °F, bringing to room temperature before using (23 °C/73 °F).
- Do not expose restorative materials to elevated temperatures or intense light.
- If using unit-dose tips, which are intended for single patient use, discard after each patient. If using syringes, always replace cap after each use.

\* Vita is a registered trademark of VITA ZAHNFABRIK H. Rauter GmbH & Co. KG Germany

	LICHTHÄRTEND MIKRO-HYBRID KOMPOSIT
	FOTOCURADO COMPUESTO MICRO HÍBRIDO
	PHOTOPOLYMÉRISATION COMPOSITE MICRO-HYBRIDE
	FOTOPOLIMERIZZAZIONE COMPOSITO MICROIBRIDO
	KIETINIMAS ŠVIESA MIKROHIBRIDINIS KOMPOZITAS
	UTWARDZANY ŚWIATŁEM KOMPOZYT MIKROHYBRYDOWY
	FOTOPOLIMERIZABIL MATERIAL COMPOZIT MICROHIBRID
	IŞIKLI KÜR MIKROHİBRİD KOMPOZİT
	PM-3913 Rev010

DENTAL  
TECHNOLOGIES  
6901 N. HAMLIN AVENUE  
LINCOLNWOOD, ILLINOIS 60712 USA  
www.dentaltech.com

CE  
0459

EMERGO EUROPE  
Prinsesgracht 20  
2514 AP, The Hague  
The Netherlands