ES

Descripción:

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EN

Description:

Alpha III™ NANO is a light cure, fluoride containing, Bis-GMA resin-based nano-hybrid composite containing fluorescence to provide duplication of the natural tooth surface. **Alpha III™ NANO** consists of 62.63% by volume of specifically selected inorganic fillers which range in particle size from 0.18 micron to 2 micron that provide exceptional wear resistance and radiant polishability. **Alpha III™ NANO** is available in 6 Vita® shades (A1, A2, A3, A3.5, B1, and B2). All shades are radio-opaque for easy identification. **Alpha III™ NANO** is packaged in a syringe, or single-dose capsules.

Indications for Use:

Alpha III™ NANO is recommended for use in all anterior and posterior restorations, including occlusal surfaces.

Contraindications:

Alpha III™ NANO is contraindicated in patients with a history of severe allergic reaction to methacrylate resins.

Warnings:

- Avoid contact of etchant with soft tissue, eyes, and skin. **Alpha-Etch37®** contains phosphoric acid, which is corrosive. In case of contact with soft tissue, eyes or skin immediately flush with water for at least 15 minutes. Seek medical attention.
- Do not use in the presence of eugenol based materials as complete setting may not occur.

General Precautions:

- U. S. Federal Law restricts this device to sale by or on the order of a dentist.
- Review instructions and SDS (www.dentaltech.com/resource) prior to beginning treatment. Use only as directed.
- If the unit dose packaging was purchased, the unit dose container is disposable and is intended for single patient use only.
- Do not expose to elevated temperatures or intense light.

Adverse Reactions:

Allergic or pre-existing skin, eye and respiratory conditions may be aggravated by exposure to this product.

Instructions for use:

- Prepare the tooth. Isolate and dry the tooth thoroughly with oil-free compressed air.
- If using **Alpha-Etch37®**, apply etchant to the prepared enamel for 15 seconds, or if applied to dentin etch for 10 seconds. Lightly agitating the etchant during its application can enhance its action. Never leave etchant on tooth for more than 60 seconds.
- Rinse thoroughly with water.
- If using a bonding adhesive, follow the manufacturer's instructions regarding indications for use, warnings and general precautions, and instructions for use.
- Select the appropriate shade of **Alpha III™ NANO**. Apply **Alpha III™ NANO** directly onto the prepared tooth in increments no larger than 2mm thickness.
- Using a visible light-curing unit (minimum output of 600 mW/cm2) cure each incremental layer for at least 40 seconds.
- Finish and polish the restoration using standard techniques. A composite surface sealant can also be used to improve surface smoothness and seal micro cracks and surface imperfections caused by the finishing procedures.

Storage:

- Recommended storage conditions are between 2°C/36°F to 27°C/80°F, bringing to room temperature before using (23°C/73°F). Do not freeze. See outer package for expiry date.
- Do not expose to elevated temperatures or intense light.
- Keep container closed when not in use.

* Vita is a registered trademark of VITA ZAHNFABRIK H. Rauter GmbH & Co. KG Germany

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