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Description:

Alpha-Crown® is a self-cure polymer-based temporary crown & bridge material with fluorescence which allows for naturally aesthetic, life-like temporaries. The **Automix** packaging design provides automatic mixing from the dual chambered cartridge to greatly reduce in-office procedure time. **Alpha-Crown®** contains 23% by volume of a blend (barium-alumino-silicate, quartz and silica) of radiopaque fillers, with a particle size range of 1 to 20 micron, for easy identification in radiographs. **Alpha-Crown®** is available in four Vita® shades (A1, A2, A3, B1). Less than 6°C generated during polymerization from a 100mm³ sample of **Alpha-Crown®**.

Indications for Use:

- Alpha-Crown®** is intended for use in temporary crowns and bridges, temporary inlays and onlays, and long term temporaries.
- Fulfills the requirements of ISO 10477:2004
- ISO 10477 does not cover the application of these polymer-based crown and bridge materials for stress-bearing areas in posterior teeth. Self-Cure Temporary Crown & Bridge Material is a Type 1, Class 1 polymer-based crown & bridge material.

Contraindication:

Alpha-Crown® is contraindicated in patients with a history of severe allergic reaction to acrylic resins.

General Precautions:

- U.S. Federal Law restricts this device to sale by or on the order of a dentist.
- Review instructions and SDS (www.dentaltech.com/resource) prior to beginning treatment. Use only as directed.
- When using a new cartridge for the first time, extrude a "pea" size amount of material to ensure an even flow and discard. Be sure to avoid cross contamination of the pastes as setting within the cartridge may occur.
- Mixing tips are disposable and intended for single patient use.
- Always leave mixing tip in place after use. It will serve as the cartridge cap.

Adverse Reaction:

Alpha-Crown® contains polymerizable monomers which can cause skin sensitization (allergic contact dermatitis) in susceptible individuals. Wash thoroughly with soap and water after contact. If such reaction occurs, discontinue use.

Instructions for Use:

- Prior to tooth preparation or extraction, take a primary impression using silicone, alginate or thermoplastic impression material. To provide bulk for later finishing of the temporary restoration, trim the interdental gingival excess from the impression. In molar areas with absent teeth, the carving in the interproximal areas should create a ridge for pontic contact between the prepared teeth.
- Perform any necessary extractions and prepare abutting teeth.
- Dry the prepared abutments. If desired, the prepared abutments and tissue can be lightly coated with Petroleum Jelly for easier separation.
- Air-dry the impression and begin dispensing **Alpha-Crown®** into the impression or matrix filling from the bottom (occlusal surface) to the top (gingival). Keeping the mixing tip immersed in the material while dispensing to prevent air bubbles.
- Insert the impression or matrix, filled with **Alpha-Crown®** within 45 seconds over the prepared teeth using steady, even pressure and hold firmly in place.
- Remove from the mouth after approximately 1 minute and 45 seconds to 2 minutes and 15 seconds (from the start of mixing) has elapsed. At this stage the material will be elastic and can be easily removed without deformation.
- Alpha-Crown®** will set in 3 additional minutes. For a quicker setting time, material can be placed within a cup of hot water (50°C/122°F) for 30 seconds.
- Remove the material from the impression or matrix. Removal of the oxygen inhibition smear layer can be easily removed by wiping with alcohol.
- If trimming is necessary, the temporary materials can be shaped with slow acrylic burs and disks. Adjust occlusion and embrasures in the mouth using high-speed finishing diamond burs. Polish with conventional materials.
- Cement the finished temporary restoration with **Alpha-TCEM™ NE** temporary cement.

Performing Repairs to Temporary Restorations:

Voids or adding to short margins:

- Rinse and air dry the temporary restoration.
- Roughen the surface to be repaired with a diamond burr.
- Remove all surface dust and wipe with alcohol.
- Apply bonding adhesive according to manufacturer's instructions.
- Apply matching shade of **Alpha-Flow®** light cure flowable composite onto repair and light cure according to the instructions.
- Finish and adjust as necessary.



Fracture:

- Repeat steps 1 -4 above.
- Dispense a fresh mix of **Alpha-Crown®** onto all surfaces.
- Hold together for 2 minutes and 30 seconds.
- Wait 5 minutes and 30 seconds before trimming or removing excess material.

Storage:

- Refrigeration is recommended. Recommended storage conditions are between 2°C/36°F to 27°C/80°F, bring to room temperature before using (23°C/73°F). Do not freeze. See outer package for expiry date.
- Keep container closed when not in use.

* Vita is a register trademark of VITA AZAHNFABRIK H. Rauter GmbH & Co. KG Germany.

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- „Vita“ yra registruotasis bendrovės „VITA AZAHNFABRIK H. Rauter GmbH & Co. KG“, Vokietija, prekės ženklas.