



# ALPHA-EASE™

## DESENSITIZER


ALPHA-EASE™ PROVIDES ↓

- FLUORIDE
- ANTIMICROBIAL
- ONE STEP PROCEDURE
- COMPATIBLE WITH ANY RESTORATIVE MATERIAL
- 91% CLINICAL RATING



### Desensitizer

ALPHA-EASE™ is a Hema/Glutaraldehyde desensitizer indicated for use for treating and preventing dentinal hypersensitivity. Simple and easy to use, can be applied to moist surface prior to the placement of restorative material. It provides a physiological seal within the dentinal tubules by polymerization with the proteins in the tubules to prevent leakage. The material provides virtually immediate relief from pain caused by hypersensitivity and is formulated with an antimicrobial to inhibit the growth of bacteria.



### DENTAL TECHNOLOGIES, INC.

6901 N. HAMLIN AVENUE  
LINCOLNWOOD, ILLINOIS 60712 USA

1.800.835.0885 TOLL FREE [US/CANADA]  
1.847.677.5500 PHONE  
1.847.677.5502 FAX