



ALPHA-CORE™ DC

DUAL CURE CORE BUILD UP MATERIAL

ALPHA-CORE™ DC PROVIDES ↓

- COST EFFECTIVENESS
- CUTS LIKE DENTIN
- RADIOPAQUE
- NON-SLUMPING
- CONTAINS FLUORIDE



Dual Cure Core Build Up Material

ALPHA-CORE™ DC is a highly filled radiopaque, dual cure resin based composite core build up material in an automix cartridge. The material is designed specifically for core build-ups and cementation of posts and pins. The material exhibits ideal mechanical properties with high compressive and tensile strength yet maintaining a dentin like feel when cutting. Its ideal handling characteristics allows for direct placement of the material with Intra Oral tips.

Alpha-Core™ DC is available in two shades: a Natural Shade (A2) for use under ceramic crowns and Blue Shade to distinguish core material from the tooth structure.

Compressive Strength	296 MPa
Flexural Strength	143 MPa
Diametral Tensile Strength	39 MPa
Depth of Cure/Light Cure @40 seconds	4mm
Depth of Cure/Light Cure @20 seconds	2mm

<p>ALPHA-CORE™ DC</p> <p>REF122.0003.002 NATURAL KIT</p> <p>REF122.0004.002 BLUE KIT</p> <p>KIT CONTAINS:</p> <ul style="list-style-type: none"> ■ 50g CARTRIDGE ■ 10 AUTOMIX TIPS ■ 10 INTRA ORAL TIPS
--

DENTAL TECHNOLOGIES, INC.

6901 N. HAMLIN AVENUE
LINCOLNWOOD, ILLINOIS 60712 USA

1.800.835.0885 TOLL FREE [US/CANADA]
1.847.677.5500 PHONE
1.847.677.5502 FAX

