



# ALPHA-CORE® SC

## SELF CURE RESIN BASED COMPOSITE FOR CORE BUILD UPS

ALPHA-CORE® SC PROVIDES ↓

- CUTS LIKE DENTIN
- EXCELLENT MECHANICAL PROPERTIES
- RADIOPAQUE
- NON-SLUMPING FORMULA
- CONTAINS FLUORIDE



### Self Cure Resin Based Composite For Core Build Ups

Alpha-Core® SC is a self cure composite specifically designed for core build ups and cementations of posts and pins. The material exhibits ideal mechanical properties with high compressive and tensile strength yet maintaining a dentin like feel when cutting. its ideal handling characteristics allows the materials to be placed with a syringe and depending on the shade of the base paste it can be used to match or contrast the color.



Compressive Strength	255 MPa
Flexural Strength	110 MPa
Diametral Tensile Strength	50 MPa

<p><b>ALPHA-CORE® SC</b></p> <p>REF 120.0025.002 WHITEKIT  REF 120.0026.002 NATURALKIT  REF 120.0027.002 BLUEKIT</p> <p>KIT CONTAINS:</p> <ul style="list-style-type: none"> <li>■ 14g CATALYST</li> <li>■ 14g BASE</li> <li>■ SPATULAS</li> <li>■ MIXING PADS</li> </ul>
---

### DENTAL TECHNOLOGIES, INC.

6901 N. HAMLIN AVENUE  
LINCOLNWOOD, ILLINOIS 60712 USA

1.800.835.0885 TOLL FREE [US/CANADA]  
1.847.677.5500 PHONE  
1.847.677.5502 FAX

