



ALPHA-GLAZE®

LIGHT CURE SURFACE SEALANT

ALPHA-GLAZE® PROVIDES ↓

- WEAR RESISTANT FINISH
- CLEAR IN COLOR
- COMPATIBLE W/ MOST RESIN BASED COMPOSITES
- IDEAL VISCOSITY
- CONTAINS FLUORIDE



Light Cure Surface Sealant

Alpha-Glaze® Surface Sealant is a fluoride containing light cure unfilled resin with low viscosity to be used as a finishing sealant for composites and resin cements. It is formulated to penetrate and seal the micro cracks and surface imperfections that are created by the finishing procedures. It can also be used to treat the enamel surface prior to applying Orthodontic Appliances.



ALPHA-GLAZE®

REF403.0101.002 ALPHA-GLAZE®

- 1 X 7mL Surface Sealant with Fluoride

DENTAL TECHNOLOGIES, INC.

6901 N. HAMLIN AVENUE
LINCOLNWOOD, ILLINOIS 60712 USA

1.800.835.0885 TOLL FREE [US/CANADA]
1.847.677.5500 PHONE
1.847.677.5502 FAX

