



# ALPHA-PRO®

## PROPHYLAXIS PASTE

ALPHA-PRO® PROVIDES ↓

- EXCELLENT STAIN REMOVAL
- SPLATTER-FREE FORMULA
- COLOR CODED TO DISTINGUISH GRIT
- GREAT FLAVORS
- CHOICE OF UNIT-DOSE CUPS OR JARS



### Prophylaxis Paste

ALPHA-PRO® is a blend of cleaning and polishing agents containing 1.23% Fluoride Ion. This splatter-free formula rinses easily and completely. Medium Grit is recommended for normal cleaning and polishing, and plaque removal. Coarse Grit is recommended for extra heavy stains, including tobacco and heavy plaque deposits. ALPHA-PRO® is available in two great flavors, Mint and Cherry. It is available in 100 g jar, 340 g jar and or 2 g unit-dose cups (200 cups per jar). ALPHA-PRO® is also available with no Flavor and no Fluoride in Medium Grit (Ref. 560.1301.002) which can be used before bonding procedures to improve bond strength.



#### ALPHA-PRO®

NET WT. 400 GRAMS  
(200 SINGLE USE CUPS)

- MEDIUM No Fluoride, No Flavor/REF560.1301.002
- MEDIUM MINT/REF560.0204.002
- MEDIUM CHERRY/REF560.0304.002
- COARSE MINT/REF560.1204.002
- COARSE CHERRY/REF560.0604.002

NET WT. 340 GRAMS  
(12 oz. JAR)

- MEDIUM MINT/REF560.0202.002
- MEDIUM CHERRY/REF560.0302.002
- COARSE MINT/REF560.1202.002
- COARSE CHERRY/REF560.0602.002

NET WT. 100 GRAMS  
(100 g JAR)

- MEDIUM MINT/REF560.0203.002
- MEDIUM CHERRY/REF560.0303.002
- COARSE MINT/REF560.1203.002
- COARSE CHERRY/REF560.0603.002

#### DENTAL TECHNOLOGIES, INC.

6901 N. HAMLIN AVENUE  
LINCOLNWOOD, ILLINOIS 60712 USA

1.800.835.0885 TOLL FREE [US/CANADA]  
1.847.677.5500 PHONE  
1.847.677.5502 FAX