



ALPHA-DENT® LIGHT CURE COMPOSITE

LIGHT CURE HYBRID COMPOSITE

ALPHA-DENT® PROVIDES ↓

COST EFFECTIVENESS

CONTAINS FLUORESCENCE

RADIOPAQUE

NON-SLUMPING

EIGHTEEN SHADES



LIGHT-CURE HYBRID COMPOSITE

Alpha-dent® Light Cure Composite restorative is radiopaque. It has the versatility of the Hybrid composite while measuring high compressive strength with excellent polishability. This composite is appropriate for almost every type of restoration. The material is available with the Exact Match System, which uses a specially designed computer program to match the composite color and translucency to the VITA®¹ shade guide.

Compressive Strength	225 MPa
Flexural Strength	95 MPa
Diametral Tensile Strength	47 Mpa

Light Cure Composite w/ Liquid Etchant System

Kit Contains: 7 x 4.5g Composite Syringes . 3ml Light Cure Bonding Resin . 7.5ml Etchant Liquid . Brushes . Brush Holder . Dispensing Pad

REF Code Product

- 101.0000.001 Hybrid Composite Standard Kit+
- 101.0001.001 Hybrid Composite Custom Kit++

Light Cure Composite w/ Etchant Gel System

Kit Contains: 7 x 4.5g Composite Syringes . 7ml Light Cure Bonding Resin . 2 x 4g Etchant Gel Syringes & Tips . Brushes . Brush Holder . Dispensing Pad

REF Code Product

- 101.0002.001 Hybrid Composite Standard Kit+
- 101.0003.001 Hybrid Composite Custom Kit++

Light Cure Composite 2-Syringe Kit w/ Liquid Etchant

Kit Contains: 2 x 4.5g Composite Syringes . 3ml Light Cure Bonding Resin . 3ml Etchant Liquid . Brushes . Dispensing Pad

REF Code Product

- 101.0008.001 Hybrid Composite Custom Kit++

Light Cure Composite 7-Syringe Kit (no accessories)

Kit Contains: 7 x 4.5g Composite Syringes

REF Code Product

- 101.0006.001 Hybrid Composite Standard Kit+
- 101.0007.001 Hybrid Composite Custom Kit++

Light Cure Composite in Unit Dose w/ Bonding

Kit Contains: 7 sets x 4.5g (15 x .30g) Composite Unit-Dose . 2 x 3ml Light Cure Bonding Resin . 2 x 4g Etchant Gel Syringes & Tips . Brushes . Brush Holder . Dispensing Pad

REF Code Product

- 101.0004.001 Hybrid Composite Standard Kit+
- 101.0005.001 Hybrid Composite Custom Kit++

Light Cure Composite 4-Syringe Kit w/ Liquid Etchant

Kit Contains: 4 x 4.5g Composite Syringes . 3ml Light Cure Bonding Resin . 7.5ml Etchant Liquid . Brushes . Dispensing Pad

REF Code Product

- 101.0010.001 Hybrid Composite Custom Kit++

+ Shades A1, A2, A3, B2, B3, C3, D4
++ Choice of shades

DENTAL TECHNOLOGIES, INC.

6901 N. HAMLIN AVENUE
LINCOLNWOOD, ILLINOIS 60712 USA

1.800.835.0885 TOLL FREE [US/CANADA]
1.847.677.5500 PHONE
1.847.677.5502 FAX