



**DENTAL  
TECHNOLOGIES**

www.dentaltech.com

# ALPHA-BOND®

ACETONE BASED BONDING ADHESIVE

ALPHA-BOND® PROVIDES ↓


SINGLE COMPONENT LC ADHESIVE

5TH GENERATION FORMULA

DIRECT & INDIRECT BONDING

DUAL CURE ACTIVATOR AVAILABLE

GREAT BOND STRENGTH



## Acetone Based Bonding Adhesive

Alpha-Bond® is a Single Component Acetone Based Adhesive with a Dual Cure Activator (purchase separately) formulated to bond to multiple substrates for both direct and indirect bonding procedures. It contains the hydrophilic monomer PMGDM that allows penetration into the tubules, intertubules and peritubular dentin. For a single step light cure adhesive, use the bottle marked adhesive. By combining the dual cure activator with the adhesive, you can confidently bond to self cure composites and areas inaccessible to light.



ALPHA-BOND®

REF 211.0001.002 6mL BOTTLE

REF 211.0201.002 3mL ACTIVATOR

## DENTAL TECHNOLOGIES, INC.

6901 N. HAMLIN AVENUE  
LINCOLNWOOD, ILLINOIS 60712 USA

1.800.835.0885 TOLL FREE [US/CANADA]

1.847.677.5500 PHONE

1.847.677.5502 FAX