




ALPHA-DENT® COPAL CAVITY VARNISH

ALPHA-DENT® PROVIDES ↓

- INTERMEDIARY VARNISH
- FORMULATED WITHOUT CHLOROFORM
- PROTECT EXPOSED DENTIN
- REDUCE FLUID LEAKAGE
- 2 YEAR SHELF LIFE



Cavity Varnish

Alpha-Dent® Cavity Varnish & Solvent is formulated without chloroform, unlike other products on the market. It is used to form a barrier to protect exposed dentin and reduce fluid leakage around the tooth restoration interface.



ALPHA-DENT®

REF 407.0000.001 COPAL CAVITY KIT

KIT CONTAINS:

- 1 X 14mL VARNISH BOTTLE
- 1 X 14mL SOLVENT BOTTLE

DENTAL TECHNOLOGIES, INC.

6901 N. HAMLIN AVENUE
LINCOLNWOOD, ILLINOIS 60712 USA

1.800.835.0885 TOLL FREE [US/CANADA]
1.847.677.5500 PHONE
1.847.677.5502 FAX