




# ALPHA-BOND® SE

## SELF ETCH ADHESIVE

ALPHA-BOND® SE PROVIDES ↓

- ELIMINATES THE ETCHING STEP
- 7TH GENERATION FORMULA
- DIRECT & INDIRECT BONDING
- GREAT BOND STRENGTH
- MINIMIZES POST-OPERATIVE SENSITIVITY



### 7th Generation Self Etch Adhesive

Alpha-Bond® SE is a self-etch adhesive formulated to bond to multiple substrates for both direct and indirect bonding procedures. It's a universal adhesive that has the bond strength and durability of a total-etch adhesive while eliminating the etching step.



### TECHNICAL DATA

Bond Strength to Dentin/Enamel 37 MPa

BRAND	SHEAR BOND STRENGTH TO DENTIN/ENAMEL*
ALPHA-BOND® SE	37 MPa
COMPETING BRAND 1	23 MPa
COMPETING BRAND 2	25 MPa
COMPETING BRAND 3	31 MPa

\*Testing completed by The Dental Advisor. Go to www.dentaladvisor.com for full report.



### DENTAL TECHNOLOGIES, INC.

6901 N. HAMLIN AVENUE  
LINCOLNWOOD, ILLINOIS 60712 USA

1.800.835.0885 TOLL FREE [US/CANADA]  
1.847.677.5500 PHONE  
1.847.677.5502 FAX