



DENTAL TECHNOLOGIES


www.dentaltech.com

ALPHA-BOND®

ACETONE BASED BONDING ADHESIVE


ALPHA-BOND® PROVIDES ↓

- SINGLE COMPONENT LC ADHESIVE
- 5TH GENERATION FORMULA
- DIRECT & INDIRECT BONDING
- DUAL CURE ACTIVATOR AVAILABLE
- GREAT BOND STRENGTH



Acetone Based Bonding Adhesive

Alpha-Bond® is a Single Component Acetone Based Adhesive with a Dual Cure Activator (purchase separately) formulated to bond to multiple substrates for both direct and indirect bonding procedures. It contains the hydrophilic monomer PMGDM that allows penetration into the tubules, intertubules and peritubular dentin. For a single step light cure adhesive, use the bottle marked adhesive. By combining the dual cure activator with the adhesive, you can confidently bond to self cure composites and areas inaccessible to light.



DENTAL TECHNOLOGIES, INC.

6901 N. HAMLIN AVENUE
LINCOLNWOOD, ILLINOIS 60712 USA

1.800.835.0885 TOLL FREE [US/CANADA]
1.847.677.5500 PHONE
1.847.677.5502 FAX