




# ALPHA III™ NANO

## LIGHT CURE NANO-HYBRID COMPOSITE

### PROVIDES ↓

- BRILLIANT AESTHETICS
- REDUCED POLYMERIZATION SHRINKAGE
- EFFORTLESS HANDLING
- IMPROVED RADIOPACITY

Alpha III™ NANO is a light cure, Bis-GMA resin-based nano-hybrid composite containing fluorescence to provide duplication of the natural tooth surface. Alpha III™ NANO consists of 62.63% by volume of specifically synthesized inorganic fillers which range in particle size from 0.18 micron to 2 micron that provide exceptional wear resistance and radiant polishability.



### AVAILABLE SHADES

- AI, A2, A3, A3.5, B1, B2

### DENTAL TECHNOLOGIES, INC.

6901 N. HAMLIN AVENUE  
LINCOLNWOOD, ILLINOIS 60712 USA

1.800.835.0885 TOLL FREE [US/CANADA]  
1.847.677.5500 PHONE  
1.847.677.5502 FAX