

TECHNICAL DATA SHEET

Product:	Alpha-Dent® Tray Adhesive
Product Description:	Alpha-Dent® Tray Adhesive is used to obtain the most accurate impressions possible by adhering impression materials to tray. Since impression materials do not have adhesive properties, it is advisable to use Alpha-Dent® Tray Adhesive on all trays before taking an impression.
Indications for Use:	Alpha-Dent® Tray Adhesive is used for obtaining accurate impressions by adhering impression materials to trays.
Appearance:	Blue viscous liquid.
Curing Mechanism:	Not applicable.
Shelf-Life:	36 months
Filler Type:	None
Main Composition (%w/w):	
	Isopropyl Alcohol 30-50
	Polyamide Resin 50 – 70
	Colorant < 0.1%
Physical Properties: (Typical Values)	
	Appearance viscous, slightly turbid liquid
	Color blue
	Specific Gravity 0.86
	Refractive index 1.438