

TECHNICAL DATA SHEET

- Product:** Alpha-Dent[®] Orange Solvent
- Product Description:** Alpha-Dent[®] Orange Solvent is a mild orange-based solvent used for the cleaning and removal of cements and pastes from the patient's skin and tissue, as well as cleaning mixing slabs and instruments.
- Indications for Use:** Designed to be used for cleaning and removal of cements and pastes from the patient's skin and tissue, as well as for cleaning mixing slabs and instruments.
- Appearance:** Clear, light yellow liquid.
- Curing Mechanism:** Not applicable.
- Shelf-Life:** 48 months
- Filler Type:** None
- Main Composition (%w/w):**

Natural and artificial oils	100
-----------------------------	-----

Physical Properties: (Typical Values)

Appearance	clear liquid
Color	light yellow
Odor	lemon
Specific Gravity	0.84
Refractive index	1.467

Regulatory Classification:

This product is in conformity with the following standards or other normative documents:

Code of Federal Regulations 21 CFR §701.3(a) – Cosmetic Labeling