

TECHNICAL DATA SHEET

Product:	Alpha-Dent® Alginate Flavorings	
Product Description:	Highly concentrated liquid flavorings that add a unique and pleasant flavor to alginate impression materials.	
Indications for Use:	Adding flavor to alginate impression materials.	
Appearance:	Clear colorless to dark-colored liquid.	
Curing Mechanism:	Not applicable.	
Shelf-Life:	30 months	
Filler Type:	None	
Main Composition (%w/w):		
	Solvent	80-90
	Sodium Saccharin	<5
	Surfactant	<1
	Flavorings:	10 – 20
	- Propylene Glycol	
	- Alcohol	
	- Water	
	- N&A Flavors	
Physical /Mechanical Properties:		
	Appearance:	Medium viscosity liquid
	Color:	Characteristic of the given flavor
	Flavor:	Characteristic of the given flavor
	Specific gravity:	0.99 to 1.09