

# TECHNICAL DATA SHEET

<b>Product:</b>	Alpha-Dent® Porcelain Conditioner	
<b>Product Description:</b>	Alpha-Dent® Porcelain Conditioner is a silane coupling agent used for increasing the bond strength between a ceramic appliance and enamel.	
<b>Indications for Use:</b>	Alpha-Dent® Porcelain Conditioner is used to increase the bond strength between ceramic appliance and enamel.	
<b>Appearance:</b>	Clear colorless liquid.	
<b>Curing Mechanism:</b>	Not applicable	
<b>Shelf-Life:</b>	30 months	
<b>Filler Type:</b>	None	
<b>Main Composition (%w/w):</b>	Components:	% w/w:
	Isopropyl Alcohol	> 80
	Silane	< 20
<b>Physical /Mechanical Properties:</b>	Refractive Index	1.376 to 1.380
	Specific Gravity	0.760 to 0.800