

TECHNICAL DATA SHEET

Product: Alpha-Dent® Orange Solvent

Product Description: Alpha-Dent® Orange Solvent is a mild orange-based solvent used for the cleaning and removal of cements and pastes from the patient's skin and tissue, as well as cleaning mixing slabs and instruments.

Indications for Use: Designed to be used for cleaning and removal of cements and pastes from the patient's skin and tissue, as well as for cleaning mixing slabs and instruments.

Appearance: Clear, light yellow liquid.

Curing Mechanism: Not applicable.

Shelf-Life: 48 months

Filler Type: None

Main Composition (%w/w):

Natural and artificial oils	100
-----------------------------	-----

Physical Properties: (Typical Values)

Appearance	clear liquid
Color	light yellow
Odor	lemon
Specific Gravity	0.84
Refractive index	1.467