

TECHNICAL DATA SHEET

Product:	Alpha-Pro® White Varnish
Product Description:	Alpha-Pro® White Varnish contains 5% sodium fluoride. It is used in the treatment of dentinal and post operative sensitivity. Alpha-Pro® White Varnish, when applied, practically becomes invisible, unlike traditional varnishes that cause unsightly yellow coloring of the teeth. Alpha-Pro® White Varnish sets immediately on contact with saliva, leaving a smooth, thin layer. It is available in bubblegum, melon, mint, and strawberry flavors.
Indications for Use:	For treatment of dentinal and post operative sensitivity.
Appearance:	White to off-white, medium viscosity liquid.
Curing Mechanism:	N/A
Shelf-Life:	30 months
Filler Type:	N/A
Main Composition (%w/w):	
	Colophony: 45 - 60
	Ethyl Alcohol: 20 - 35
	Sodium Fluoride: 5
	Emulsifier: 2 - 7
	Flavor: 2 - 5
	Xylitol: <1
	Sweetener: <1
Physical /Mechanical Properties:	
	Appearance: Medium viscosity liquid
	Color: White to off-white
	Odor: Characteristic of ethanol and specific flavor
	Flavor: Bubblegum, melon, mint or strawberry
	% Sodium Fluoride: 5