

TECHNICAL DATA SHEET

- Product:** Alpha-Dent® Light Cure Band Adhesive
- Product Description:** Alpha-dent® Light Cure Band Adhesive is a single paste resin-based, fluoride containing band adhesive. The single paste adhesive requires no mixing; resulting in reducing adhesive waste, consistent adhesive performance and no working time constraints. Alpha-dent® Light Cure Band Adhesive's distinctive blue color allows you to easily see the adhesive during banding, making clean up easy. When used with a high performance curing light, cure time is 20 seconds.
- Indications for Use:** Alpha-dent® Light Cure Band Adhesive is a light-cured, resin-based, fluoride releasing paste used for cementation of orthodontic appliances to teeth.
- Appearance:** Blue, opaque, medium viscosity paste.
- Curing Mechanism:** Light cured with dental curing light with a minimum intensity of 600 mW/cm² in the 400-500 nm range of the visible light spectrum.
- Shelf-Life:** 30 months at 25°C and 60% relative humidity.
- Filler Type:** Fluoride containing barium-alumino-silicate glass.
- Filler by Volume: 47%
Filler Particle Range: 0.34 – 2.19 micron
Average Filler Particle Size: 0.77 micron

Main Composition (%w/w):

Inorganic filler	60 – 70
Methacrylate monomers	25 - 35
Thickeners	2 - 5
Accelerator	<1
Photoinitiator	<1
UV Stabilizer	<1
Colorant	<1

Physical /Mechanical Properties: (Typical Values)

Consistency	0.50 mm
Sensitivity to ambient light	>60sec.
Depth of cure [actual]	5.11 mm/20sec.
Bond strength to enamel	5.50 MPa

DENTAL TECHNOLOGIES, INC.

6901 N. HAMLIN AVENUE
LINCOLNWOOD, ILLINOIS 60712 USA

1.800.835.0885 TOLL FREE [US/CANADA]
1.847.677.5500 PHONE
1.847.677.5502 FAX